
CATALOGUE 2025

Arabia for science & Technology

PHYSICS EXPERIMENTS

25 YEARS OF EXPERIENCE

ARABIA4SCIENCE.COM

01016322246 - 01117293405
01099973559 - 01065101051
0244050900 - 0244050919

info@arabia4science.com
sales@arabia4science.com

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Aim of the experiment:

Practice on determination of different dimensions using Vernier caliper, micrometer and spherometer.

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Aim of the experiment:

- The objective of simple pendulum experiment is to Study the motion of simple pendulum, to determine the acceleration due to gravity 'g'.

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
Aim of the experiment:

- Determining the force constant or elastic constant of a spring.

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
Aim of the experiment:

Determination of the acceleration due to gravity by means of a compound pendulum.

Harmonic motion of the compound pendulum is no simpler like the simple harmonic one.

The motion of the compound pendulum is studied here with the wall mounted metallic rod with the ability to change the two arms lengths of the rod According to the relation

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Aim of the experiment:

- Apparatus for measuring the time it takes for a ball to fall a certain distance using a digital timer. Very easy to set up and use but nevertheless highly accurate. Includes 3 steel balls.
- To determine gravitational acceleration by studying the velocity of a falling object as a function of time

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OPTIONAL

PHOTOGATES CAN BE CHANGED TO ELECTRICAL RELAYS

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Aim of the experiment:

- Study the projectile motion.
- Determine the maximum height (h) and the maximum range (S) of projection as a function of the angle of inclination (φ).

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
Aim of the experiment:

The aim of that experiment is determination of the moment of inertia of a flywheel.

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
- Determining the moment of inertia, I of the disc
- Demonstrate precession and nutation.

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Component
Disc
Gyroscope axle
Bobbin
Set of accessories U52010
Tripod stand
Stand rod
Counterweight, large
Counterweight, small
Adjusting screw
Scale
Spirit level

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- Maxwell wheel with large moment of inertia for demonstrating the conversion of energy in the conversion of kinetic energy to potential energy and vice versa. Includes suspending rod and adjustable suspension mechanism.

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
Aim of the experiment:

- The purpose of this experiment is to use the force table to experimentally determine the force which balances three other forces.

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Aim of the experiment:

The law of the lever follows from the equilibrium of moments, which works for all three classes of lever.

It represents the physical basis for all kinds of mechanical transmission of force.

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**Aim of the experiment:**

- An experimental demonstration of Newton's second law.

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**Aim of the experiment:**

- This experimental setup is used to study many of the mechanics laws like: uniform motion laws, average and instantaneous speed, uniform accelerated motion, elastic collision, frictionless motion, and many more.

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
Aim of the experiment:

- To measure the coefficient of static friction.
- Determination the downward force as a function of the angle of inclination.
- Comparison of static and dynamic friction.
- Measurement of how dynamic friction depends on the area in contact.
- Measurement of how dynamic friction depends on the combination of materials.
- Measurement of how dynamic friction depends on the perpendicular force between the two surfaces (normal force).

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It is a comprehensive mechanical experiment device that can perform multi-purpose experiments including levers, pulleys and weights, balance force, axels, pendulum motion, and slope experiments

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Component	
1 slope slide (200x600x100mm support and box)	3 pendulum weights
1 pulley	cart 1 piece
2 each of 2 pulleys	1 plate
20 weights of 20g	1 parallelogram beam (including 2 pulley mechanisms)
1 lever (span)	thread
1 3-stage axle wheel	some hooks

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Aim of the experiment:

- Young's modulus of a wire (for example) is a measure of the elasticity of the wire i.e. the ability of the wire to restore its initial length after removing the force applied to it.
- This experiment setup measure the young's modulus of a wire by measuring the longitudinal stress and strain of the wire according to the relation

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
- A flat bar is supported at one end. It is bent by the action of a force acting at the other end.
- The modulus of elasticity is determined from the bending and the geometric data of the bar.

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
- The rod is supported at two points. The rod is bent by a force acting at the center of the rod.
- The modulus of elasticity is determined from the bending and the geometric data of the bar.

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
Aim of the experiment:

- The torsion equipment is designed for determining the torsion coefficient and shear modulus of cylindrical metal wires by means of static measurements of the torsion angle and the force of torsion as well as dynamic measurements of periods of oscillation of a torsion pendulum.

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Pendulum disc with 4 pins attached	Cylindrical steel wire (diam.: 2 mm, length: 500 mm)
Upper cross piece	Clamping arm for cylindrical rod
Vertical rod	Fastening screws for clamping rod
Slot for upper cross piece (for cylindrical rods of length 500 mm or 300 mm)	Weights
Horizontal rod	Slot for lower cross piece
Clamping sleeve	Low crosspiece (with rubber feet)
Cylinder to fit clamping sleeve	Scaled disc
Precision dynamometer, 2 N	Precision dynamometer, 5 N

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Aim of the experiment:

To determine the acceleration due to gravity by means of a conical pendulum.

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- The instrument is used to measure the period of a pendulum as a function of the length of the pendulum and the effective component of the earth's gravitational field ("acceleration due to gravity").

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
Aim of the experiment:

- Measurement and analysis of simple harmonic rotary oscillation

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
**PROPERTIES OF MATTER AND FLUID
EXPERIMENTS**

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Aim of the experiment:

- Determination of the viscosity of a liquid.

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Aim of the experiment:

- Ostwald viscometer, also known as U-tube viscometer or capillary viscometer is a device used to measure the viscosity of the liquid with a known density.

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Aim of the experiment:

- Determine the specific weight of a solid body.
- Verification of Archimedes principal
- Determine the density of water (or any liquid sample)

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
Aim of the experiment:

- The electronic hydrostatic balance is an electronic balance specially designed for earthwork. It is equipped with a water tank, a tray, and a universal wheel. It is easy to operate.
- Used for the measurement of solid sample density and volume. It is a combination of an electronic balance and a hydrostatic device.

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- For demonstrating the hydrostatic paradox and for quantitative measurements of ground pressure. This pressure is measured through the curvature of a membrane and indicated in magnified form with the help of a lever multiplier. Compensation for comparative measurements is possible. Includes four differently shaped vessels made of glass.

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Aim of the experiment:

- simple demonstration that allows you to measure the specific gravity of an unknown oil using a u-tube manometer. It also gives students the opportunity to visualize a multi-fluid manometer.

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Aim of the experiment:

- In this experiment, we can study the lenses and mirror laws.
- We can determine the focal length of convex and concave lenses.
- We can determine the focal length of convex and concave mirrors.
- To investigate the real and virtual images formed by concave and convex lenses.

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Aim of the experiment:

- In this experiment, we can study the diffraction of light which is defined as when it passes through an edge or slit the light is diffracted.
- To determine the width of single slit and double slit by using wavelength or the opposite.

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
Aim of the experiment:

- Determination of the Wavelength of Laser Source Using Diffraction Grating

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**Aim of the experiment:**

- In this experiment, we can determine the refractive index of glass prism.
- The refractive index of glass is determined by refraction of light through the prism at minimum deviation.

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
Aim of the experiment:

- In this experiment, We can study a traveler microscope is used to determine the refractive index of glass sheets by using a classical method. Glass sheets with different thickness are used also to determine the refractive index.

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- In this experiment, we can determine the power of liquid lens.
- We can determine refractive index 'n' of a substance which is a dimensionless number that describes how light propagates through that medium.

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
- In this experiment, we can study the relationship between the sugar concentration and refractive index and see that the relation is approximately linear.
- We can determine the unknown concentration of sugar solution.
- We can see the deviation of laser beam when pass through the prism.

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Aim of the experiment:

- In this experiment, we can Measure the rotation of the plane polarized light passing through sugar solution, by knowing the sugar tube length and concentration, the specific rotation of the sugar solution can be calculated.

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Aim of the experiment:

- To observe the rotation of the plane of polarization of monochromatic light by sugar solution for determining the concentration of solution of optically active substance.

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Polari-meter Laurent's half shade
Sodium Light Source with power supply and housing
Polari-meter tube (big)
Digital balance 2dig
Glass beaker 250ml


Specification

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Circular scale: 0° - 360°
Least count: 1°
Vernier Reading: 6 min
Dimensions: Dia. 12 mm, length 200mm
Tube length: 200mm
Tube: With central bulb, metallic caps & cover glasses packed in a velvet case

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Aim of the experiment:

- In this experiment, we can see Monochromatic light falls on a mica plate perpendicular to its optic axis. At the appropriate plate thickness ($1/4$, $1/2$ plate) there is a phase shift between the ordinary and the extraordinary ray when the light emerges from the crystal.
- The polarization of the emergent light is investigated at different angles between the optic axis of the $1/4$, $1/2$ plate and the direction of polarization of the incident light.

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**Aim of the experiment:**

- In this experiment, we can verify malus' law which states that when completely plane polarized light is incident on the analyzer, the intensity I of the light transmitted by the analyzer is directly proportional to the square of the cosine of angle between the transmission axes of the analyzer and the polarizer.

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**Aim of the experiment:**

- In this experiment, we can study the interference and its type that is newton's ring is interference fringes of equal thickness.

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Aim of the experiment:

- By dividing up the wave front of a beam of light at the Fresnel bi-prism, interference is produced. The wavelength of the light source is determined from the interference patterns.

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Aim of the experiment:

- In this experiment, we can study the interference by Michelson interferometry, which is In a Michelson interferometer, light is split up into two beams by a half-silvered glass plate (amplitude splitting), reflected by two mirrors, and passed again through the glass plate to produce interference phenomena behind it.
- To determine the wavelength of the light source.

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Aim of the experiment:

- Understand this study and measure the wavelength of the helium-neon laser by applying Michelson interferometer and coherence.

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**Aim of the experiment:**

- In this experiment, we can determine the wavelength of Hg lamp.
- Most light is a mixture of different colors and the maximum for each color occurs at a different angle because each color has a different wavelength, and each element has a unique line spectrum.

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
**Aim of the experiment:**

- In this experiment, we can study the photoelectric effect, which is defined that electrons are emitted from solids, liquids, or gases when they absorb energy from light. Electrons emitted may be called photoelectrons.

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
Aim of the experiment:

- Comprehensive set of equipment with built in He-Ne laser. Top quality optical glass components allow over 47 experiments in geometric and physical optics.

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Main unit	Included a diode laser (powered by 3 AA batteries), beam expander, beam splitter, whiteboard, dial disk, etc.	
Accessories	Right angle prism	Double plane mirror
	Equilateral prism	Semi-cylindrical lens
	Concave-convex cylindrical mirror	Beam expander lens
	Diffuse reflector	Plastic optic fiber
	Plane parallel plate	Plane mirror
	Bi-convex cylindrical lens	Plano-concave cylindrical lens
	Galilean telescope	Periscope
	Plano-convex cylindrical lens	PDF manual

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Aim of the experiment:

- In this experiment, we can determine the electro optical curve of liquid crystal sample and main parameters such as threshold voltage, saturation voltage, contrast, and gradient.
- Demonstrate the basic principle of liquid crystal display.
- Some other elements are provided for polarized light experiment and verification experiment of Malus Law.

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Optical bench	Current meter
Polarizer	Power meter
Analyzer	Connecting cables
LC plate with type TN	Voltmeter
Semiconductor laser	Photo detector laser

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Aim of the experiment:

- This apparatus can demonstrate how to produce laser hologram.

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
Aim of the experiment:

- In this experiment, we can study Holographic photography which is based on the interference principle resulted from the coherent superposition of optical beams.
- to record the interference fringes between a reference beam and an object beam (reflected from an object) on a recording medium.
- The interference fringes contain the amplitude and phase information of the object beam.
- This experiment kit employs a semiconductor laser as the light source to generate reflective and transmissive holograms on the recording medium.
- a photopolymer plate. In addition, the recorded reflective hologram can be reconstructed with white light.

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Semiconductor laser.
Exposer shutter and timer.
Beam expander.
Plane mirror.
Small object.
Rad sensitive polymer plate.
Optical table.
Metallic base

ELECTRICITY AND MAGNETISM EXPERIMENTS

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**Aim of the experiment:**

- Magnetic moment affects the magnetometer by deflecting its needle. If two magnets with different magnetic moments are placed on both sides of the wooden arm, the magnetic moments of those magnets deflect the needle with a net deflection depending on the difference between the two moments.

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
**Aim of the experiment:**

- A classic method for determining the resistance is using a meter bridge method with sensitive center reading galvanometer and 1 meter wire with current balance concept between the two sides of the bridge.
- We can use this bridge to study the series and parallel connection.

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Aim of the experiment:

- In this experiment we can determine the magnetic field strength. As the current passing through a coil, magnetic field is produced and the intensity of this is proportional to the current intensity and number of coils turns. Three coils with different number of turns and adjustable power supply is used to generate the magnetic field and a compass located in the middle of the coils is used to obtain this magnetic field.

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**Aim of the experiment:**

- Recording the characteristic curve of light dependent resistance.

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Aim of the experiment:

- To demonstrate electromagnetic induction and transformer principles. Primary coil with 600 turns, secondary coil with three different turns. Coils are wound on molded pipe and are provided with 4 mm safety sockets.

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Aim of the experiment:

- In the experiment Biot-Savart law, the profile of the magnetic flux density B of current-carrying conductor loops with different radii R and the magnetic flux density of a straight wire is determined and discussed.

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
Aim of the experiment:

- In this experiment, we can study the magnetic field of an air coil in which the length L can be varied for a constant number of turns N . For the magnetic, field the relationship
- To Measure the magnetic field B of a long air coil as a function of the current I .
- To determine the magnetic field constant μ_0 .
- To measure the magnetic flux density along the axis of long coils

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Aim of the experiment:

- Helmholtz Coils are a special arrangement of coils that are placed in such a way as to generate a very uniform magnetic field between them when a current is applied. The strength of the magnetic field is proportional to the number of turns in the coils and the current applied to them. Helmholtz Coils setup consists of two field coils mounted on a track to provide a uniform magnetic field between the coils. The proper separation for Helmholtz coils is marked on the track. A pair of precision coils and 400 turns of wire are wound on heat resistant former of about 150 mm diameter.

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Aim of the experiment:

- When an electric current flow in conductor magnetic field is produced, when pass this conductor in horseshow magnet and the direction of current of wire perpendicular to horseshow magnet we can note force acting on the wire the direction of the force on wire dependent on the direction of current flow in wire this farce is called Lorentz force.

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- Apparatus is designed primarily for demonstration of the basic principles of transformers. It consists of a laminated U-core with laminated I-core, both with 30 x 25 mm cross section to form a closed core of size 105 x 140 mm. Heavy aluminum alloy stand with removable clamping rods, 2 pole pieces, each 65 x 25 x 25 mm with one end cone shaped, to accept support rod and shading ring. Poles also have a hole drilled through lengthwise to take a light beam.

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
Aim of the experiment:

- To experimentally find out the relation between the measured temperature and the thermal electromotive force generated by the real thermocouple.
- To determine Seebeck coefficient.
- To determine room temperature.

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Aim of the experiment:

- In this experiment we can Determine the dielectric constant for air ϵ_0 and ϵ for glass or acrylic sheets using two isolated parallel plates (as a capacitor) using capacitance meter.

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

Aim of the experiment:

- In this experiment we can Determine the dielectric constant for air ϵ_0 and ϵ for glass or acrylic sheets using two isolated parallel plates (as a capacitor) and high voltage power supply to charge the dielectric material.

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Aim of the experiment:

- to study the characteristics of the electric field by mapping the equipotential lines and constructing electric field lines of various charge configurations.

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Aim of the experiment:

- It's an experiment of a new process directly measuring. coulomb force by high voltage on the parallel plates electrode and an electric balance beam. It's very easy. way of verifying coulomb force

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Micrometer On Slide
Parallel Plate, $\phi 150\text{mm}$
Balance, 0~310g, 0.01g
High Voltage Connection Lead
Kilovoltage Power Supply, 15KV

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Aim of the experiment:

- The Waltenhofen's pendulum is used to demonstrate the working of an eddy current brake and its effect.

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
Aim of the experiment:

- Determine the resistance of semiconductors as a function of Temperature.
- Determine the energy gap of semiconductors.

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
Energy gap kit consists of NTC and PTC sample inside a copper or aluminium block with the thermocouple all based on one board with temperature controller with on/off switch for controlling heating

2 x Digital multi-meter.

Connection cables.

The copper cube can be changed for other materials such as (Aluminium - Stainless steel - Iron)

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Aim of the experiment:

- Learning the principle of hall effect and mapping the curve of hall element.
- Learning to measure the magnetic induction strength b using the hall effect.
- Measure the magnetic field distribution of the electromagnet coil.

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Electromagnetic excitation current: 0 ~ 500 mA is linear adjustable; adjustment fineness is <1 mA and stability is <10⁻⁵, 3-digit semi-digital voltmeter display.

Magnetic induction intensity of the magnetic field adjustable range is 0 ~ 100mT; the magnetic field is about 100mT when the current is 500 mA.

Excitation ammeter: 0 ~ 2000 mA, resolution of 1 mA error $\pm 0.5\%$;

Hall voltage measurement meter: 0 ~ 200 mV, resolution 0.1 mV error $\pm 0.5\%$; 0 ~ 2000 mV resolution 0.1 mV error $\pm 0.5\%$.

Working ammeter: range of 20 mA, resolution of 0.01 mA, error of $\pm 0.5\%$.

Gallium arsenide Hall element: the working current of 0 ~ 10 mA is linearly adjustable.

Sensitivity: 150 mV / (mA * T).

Unequal potential difference: <1 m

Two-dimensional moving ruler movement range: horizontal 0 ~ 70mm, vertical 0 ~ 30mm.

Use C-shaped electromagnet to provide magnetic field, magnetic core position can be fine-tuned, fine-tuning distance of 1mm.


The excitation current is switched by a high quality longevity 10A large current switch

The Hall element protection device is provided to avoid burning the Hall element due to the working current source and the excitation current source in the experiment, which improves the use safety of the instrument.

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**Aim of the experiment:**

- Study the hysteresis loop in ferromagnetic materials.

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THERMODYNAMICS AND GASES LAWS EXEPRIMENT

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**Aim of the experiment:**

- The Boyle-Mariotte device is used to investigate the relationship between volume and pressure of a sealed volume of air at a constant temperature. It confirms the law of Boyle-Mariotte.

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Aim of the experiment:

- Toward the end of the 1600s, the French physicist Guillaume Amontons built a thermometer since the pressure of a gas is directly proportional to its temperature at constant volume. The relationship between the pressure and the temperature of a gas is therefore known as Amontons' law.

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- This experiment used to study the physical meaning of the electrical equivalent of heat and used to determine it experimentally.
- When electric current pass-through resistance, the electric energy is converted into heat.

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Aim of the experiment:

- In this experiment we determine the specific heat capacity of solid by mixing.
- The specific heat capacity is the amount of heat required to raise the temperature of unit mass of substance 1C°.
- When two substances of different temperatures are mixed, the amount of the heat given by the first substance equals the amount of heat taken by the other one and We can use the relation:

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**Aim of the experiment:**

- Determination of the thermal conductivity of (wood, glass, and acrylic) discs by Lee's Disc method.

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Aim of the experiment:

- Determination of the thermal conductivity of metals with known physical parameters by using Searle's method.

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**Aim of the experiment:**

- Determination of the coefficient of thermal conductivity of rubber

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Aim of the experiment:

- Consider two different liquids with the same volume contained in two similar calorimeters and placed in the same atmosphere, and then the loss of heat at a second is equal for the two liquids. But if the two liquids have the same temperature and cools in the same range of temperature, they take a different time because of their difference in specific heat and the liquid with has a higher specific heat, has a lower rate of heat loss. Let the specific heat of the liquid C and for the water equal to one, then the rate of loss of heat in water is equal to the rate of loss of heat in liquid, So The quantity of heat lost by water = the quantity of heat lost by liquid.

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- In this experiment, we can study the thermal expansion of the metal.
- We can determine the thermal expansion coefficient of different metal (copper, stainless steel, and aluminum).

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$$\alpha = \frac{1}{L} \cdot \frac{\Delta L}{\Delta T} , \text{ Where: } \alpha \text{ is the expansion coefficient (the unit of } \alpha \text{ is } K^{-1})$$

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Aim of the experiment:


- In this experiment, we determine the melting point of wax from its cooling curve which is the temperature at which the substance change its state from solid to liquid. The Determination of melting point is based on latent heat which is the amount of heat required to change a unit of mass (1gm) of a substance from one state to another state without changing its temperature. So, the melting point of substance can be measured at the part where the temperature is constant with changing time.

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Aim of the experiment:

- A small amount of ice is placed in a calorimeter containing water. By knowing the masses of the ice, the water, and the calorimeter, and the resulting temperature change after the ice melts, the latent heat of fusion of ice is found.

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
**Aim of the experiment:**

- to determine the latent heat of vaporisation of water

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Aim of the experiment:

- Weld tail ends of both copper coil (Cu) and bismuth coil (Bi) or bismuth coil and antimony coil (Sb), and there occurs an electric potential difference and flows electric current when it's heated on the linked spot. It's called the Seebeck effect because of those electric current flows on a closed circuit by a voltage from a thermostatic difference.

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
Main Body
Display Equipment
Cup
Connection Lead

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**Aim of the experiment:**

- According to Stefan-Boltzmann's law, the energy emitted by a black body per unit area and unit time is proportional to the power "four" of the absolute temperature of the body. Stefan-Boltzmann's law is also valid for a so-called "grey" body whose surface shows a wavelength-independent absorption coefficient of less than one. In the experiment, the "grey" body is represented by the filament of an incandescent lamp whose energy emission is investigated as a function of the temperature.

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
Aim of the experiment:

- The Leslie's cube is used for investigating heat radiation from a hot body as a function of its temperature, color, and surface characteristics. Leslie's cube is a hollow cube with a removable lid for pouring in hot water and openings for inserting a thermometer or thermal sensor.

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Aim of the experiment:

- Peltier heat pump in symmetrical and asymmetrical operation
- Recording of temperature over time during pumping and settling phases
- Heating, cooling and electrical power, power coefficient, efficiency.
- See-beck coefficient
- Peltier effect, Joule losses, heat reflux due to the conduction of heat by Peltier element
- Heat transfer and thermal conduction coefficients
- How the thermal voltage and operating voltage depend on the temperature difference.

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Transport case	Wiring diagram sketch
Schematic view	Recess for stirrer unit
Stirrer unit	Digital thermometer for "cold side"
Connector for stirrer unit	Digital thermometer for "hot side"
Connector for power supply	Recess for stirring rods
Stirring rod "cold side"	Reservoir for "hot side"
Peltier element	Stirring rod for "hot side"
Reservoir for "cold side"	Connector for measuring current (operating current provided by Peltier element)

Aim of the experiment:

- Temperature Measuring Bench is designed to demonstrate several commonly used methods of temperature measurement and provides the means for calibration and accuracy comparisons of the different methods, so suitable measuring devices can be found for particular applications. Features are provided so that several faults commonly occurring in thermocouple systems can be demonstrated. In addition, the unit may be used to provide a temperature measuring facility for use with other experiments. The main unit comprises a bench mounting console which contains the hot water bath, furnace, ice bath, a digital resistance thermometer/thermocouple instrument, and digital thermocouple instrument and power supplies. The hot water bath comprises a Dewar flask whose lid is fitted with a combined heater stirrer unit; low water cut out to provide boiling water reference at 760 mm Hg barometric pressure. The heater unit incorporates an adjustable thermostatic temperature control and over temperature cut-out.

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
**Aim of the experiment:**

- The aim of this experiment is studying the standing waves.
- This experimental setup shows the relation between the frequency of vibration and velocity of propagation of stationary wave in a spring.

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Aim of the experiment:

The Experiment Covers the Following:

- One-Dimensional standing wave resonance of string.
- Two-Dimensional resonance using
- Chladni plates.
- Standing waves in Bohr's atomic model,
- Transverse oscillation of cantilever beam.
- Longitudinal oscillation of spiral spring

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
Aim of the experiment:

- Sound in air propagated only by longitudinal waves in which the particle's motion consist of oscillation back and forth in the direction of propagation.
- Using the properties of wave motion, the frequency, and the speed of the sound in the rod can be determined.

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Aim of the experiment:

- Resonance tube is designed to demonstrate interference of sound waves. The apparatus consists of a resonance tube with a scale which is connected via a silicone tube to another vessel with which the height of the water column can be adjusted.
- This setup used for determining velocity of sound in air

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Aim of the experiment:

- A sonometer is defined as. The device that is used for demonstrating the relationship between the frequency of the sound that is produced by the string when it is plucked and the tension, length, and mass per unit length of the string. The sound is produced in the transverse standing wave in the string.
- To study the relation between the frequency and the length of a given wire under constant tension using a sonometer.

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

**Aim of the experiment:**

- Introduction to microwaves and measure the microwave signal and penetration properties of materials.
- Reflection, diffraction, interference, and polarization of microwaves can be studied

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
Aim of the experiment:

- In the experiment two identical ultrasonic transducers are used as the transmitter and the receiver and differ only in their connection. One transducer is mounted on a measuring trolley with electric drive, while the other transducer is at rest on the lab bench. The frequency of the received signal is measured using a high-resolution digital counter.
- To determine the speed of the transducer in motion, the time Δt which the measuring trolley requires to traverse the measuring distance is measured using a stopwatch.

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
Aim of the experiment:

- Demonstrating the principle of an echo sounder.
- Determining the velocity of sound in air from the transit time of a sound pulse and the distance to the reflecting object.
- Determining distance by measuring the transit time of the sound pulse.

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Aim of the experiment:

- Sound needs a material medium with which it can enter into reciprocal action for its propagation, whereby a loss of energy occurs.
- The amplitude and so also the intensity decreases along the propagation path.
- This can be detected by measuring the signal of a receiver.

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Aim of the experiment:

- The Ripple tank provides a dramatic demonstration of the general properties of waves and propagation phenomena.
- Reflection and refraction. By using the linear dipper bar plane parallel waves can be produced. Reflection and refraction of waves can be demonstrated by using appropriate barriers in the water tank.
- The famous double slit experiment: When a plane wave encounters a barrier with two holes, these act like two-point source dippers, giving rise to the same interference pattern.
- Interference occurs when two-point source dippers generate circular waves. The distance between the sources and their frequency can be regulated.

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Ripple Tank 40 x 40 cm.	Concave Perspex Plate.
Four push Fit Legs.	Convex Perspex.
Ripple motor with eccentric cam, mounted on PVC bar.	Rectangular Perspex Plate.
Three supports for the illuminator and the motor unit.	Curved Reflector.
Motor hanging clamps.	Barrier Large.
Spring for motor support.	Barrier Small.
Lamp Unit halogen 6 Volts.	Plastic Roller Bar.

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- The periodic variations of density caused by an ultrasonic standing wave pattern in a liquid act as an optical grating for the diffraction of a monochromatic parallel light beam that is transmitted in the direction perpendicular to the direction of the ultrasound wave. From the diffraction pattern and the known wavelength of the light, it is possible to determine the sound wavelength and use that to calculate the velocity of sound in the liquid.

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Ultrasonic cw Generator
1 Test Vessel
1 Laser Diode for Debye-Sears Effect, Red
1 Laser Diode for Debye-Sears Effect, Green
1 Pocket Measuring Tape, 2 m
1 Ultrasonic Coupling Gel


Aim of the experiment:

- Large equipment set incorporating 30 student experiments for demonstrating the fundamental properties of waves using the example of 40 kHz ultrasonic waves. Stored in a tough Gragnell tray with foam inlay featuring recesses moulded to the shape of the apparatus and covered by a transparent lid. Includes CD with experiment instructions. Includes two ultrasonic transmitters, a rod-shaped microphonic sensor for recording and analysing oscillations using a standard oscilloscope and an ultrasonic pen for recording wave fronts along the desktop in the form of lines of the same phase (isophases). Many of the experiments can also be carried out without using an oscilloscope. In order to measure ultrasonic amplitudes, it is sufficient in many cases to use an analogue voltmeter for alternating current if it has a wide enough frequency range.

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1 Ultrasonic control unit	1 Concave mirrors
2 Ultrasonic transmitters, 40 kHz	2 Side pieces for double slit/reflectors
1 Ultrasonic pen	1 Centre post for double slit
1 Holder for ultrasonic pen	1 Clap for double slit
1 Holder base for ultrasonic pen	1 Ultrasonic absorber
1 Microphone probe	2 BNC cables, 1 m
2 Beam splitters	1 Cable, BNC/4-mm
3 Clamps for beam splitters	1 Plug-in power supply
1 Fresnel zone plate	

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Aim of the experiment:

- 1MHz pulse signal which is generated from the circuit board and transformed into optical signal is transmitted through a fiber-optic cable. Calculate a sonic speed with both of that pulse signal generated from the circuit board and of that optical signal transmitted through 20m fiber-optic cable.

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

**Aim of the experiment:**

- Studying the relation between the Current, voltage, and resistance.
- Verification the law of resistances in series and parallel.

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
**Aim of the experiment:**

- According to Kirchhoff's current law, the sum of currents flowing towards a point in an electrical circuit (that does not represent a capacitor plate) is equal to the sum of currents flowing away from that point.
- The algebraic sum of all the voltage around closed path is zero or in other words, the sum of the voltage drops equals to the total source voltage.

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
**Aim of the experiment:**


- The aim of the experiment is to investigate the concepts of voltage and current division.
- The sum of the voltage drops equals to the total source voltage.

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Aim of the experiment:

- Determine the Time constant of RC circuit.

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

**Aim of the experiment:**

- The objective of this experiment is to observe and measure the transient response of a series inductor-resistor, RL circuit.

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
Aim of the experiment:

- The experiment aims to study the electrical characteristics of an RLC circuit in series. Also,
- To study the relation between the input frequency f and the circuit impedance Z .

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Aim of the experiment:

- The behavior of (RC Low-Pass Filter - RC High-Pass Filter - RLC Low-Pass Filter - RLC Band-Pass and Band-Stop Filter) are studied by measuring the frequency response and step response

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
**Aim of the experiment:**

- The aim of the experiment is to compare the characteristics of red, yellow, and green light-emitting diodes.
- To estimate the wavelength of the emitted light by using the following relation.

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
Aim of the experiment:

- Study the characteristics (I-V curve) of Zener diode.

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

Aim of the experiment:

- Study the characteristics (I-V curve) of Si-&Ge diodes.

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Aim of the experiment:

- Determined Input characteristics and output characteristics and h_{fe} in CE configuration.

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The kit demonstrates the NPN transistor with various components:

Digital built-in power supplies V_{BE} (0-3 V DC) and V_{CE} (0-20 V DC)

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
Aim of the experiment:

- Determination of voltage gain
- We will learn about the working of a Transistor as a Switch. Switching and Amplification are the two areas of applications of Transistors and Transistor as a Switch is the basis for many digital circuits.

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**Aim of the experiment:**


- Recording the characteristics of a field-effect transistor.

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**Aim of the experiment:**

- To study the basic characteristics and application of the operational amplifiers

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

NUCLEAR AND MODERN PHYSICS EXPERIMENTS

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Aim of the experiment:

- To study the emission of light from a hydrogen discharge source.
- To learn the empirical formulas to characterize the pattern of spectral lines from hydrogen.
- To study and develop the Bohr Theory of the hydrogen atom.
- To measure the wavelengths of the Balmer series of visible emission lines from hydrogen.
- To learn to analyze the wavelength data to determine the Rydberg constant using the Bohr model formulation.

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Aim of the experiment:

- To observe the photoelectric effect.
- To how kinetic energy is defined using stopping potential.
- We examined each of the five component colors measure the stopping potential.
- Measuring Planck's constant.

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Aim of the experiment:

- To demonstrate the quantization of light and to measure Planck's constant (h)

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
The kit consists of a set of LEDs emitting light across a range of wavelengths (from deep blue to infrared) from wavelength of 395nm to 940nm with a built-in digital DC power supply 9V and a digital display used as voltmeter

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Aim of the experiment:

- In the Franck-Hertz experiment, an electron beam is produced by thermionic emission from a filament.
- The electrons are accelerated, pass through the vapor, and are then retarded (decelerated) by a few volts before collection at the anode.
- This all takes place in a tube contained within the Argon or Neon vapor density.

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Aim of the experiment:

The experimental instrument is an integrated experimental instrument with compact design, intuitive panel, complete functions and convenient operation

Experiments:

- Understand the general principle and usage of computer real-time measurement and control system.
- The influence of temperature, filament current and other factors on the F-H experimental curve is analyzed.
- The existence of atomic energy level is confirmed by measuring the first excitation potential of argon atoms.

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Description	Specifications
Mainbody	Display and operation with LCD screen
Power Cord	
Data Wire	
Experimental tube	Argon tube
Temperature Control Device	Control the temperature of the Argon tube

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Aim of the experiment:

- To experimentally demonstrate the concept of Millikan's oil drop experiment.
- To find the terminal velocity of the drop.
- To find the charge on a drop.

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Aim of the experiment:


- Master the method of measuring the charge of an electron using the equilibrium method.
- Demonstrate the discontinuity of the charge and determine the magnitude of the basic unit charge value.

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This compact Millikan Oil Drop Apparatus has been upgraded for increased performance and convenience and includes a permanently- installed video camera on the observing microscope, digital read-outs for the capacitor voltage and timer, and an 8” monitor for observation and control. The unit comes in a sturdy metal case, and includes silicone oil, an atomizer, and a dropper for filling the oil.

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Aim of the experiment:

- Thermionic emission of electrons
- Deflection of electrons in electron and magnetic field
- Estimate of specific charge of an electron e/m

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Aim of the experiment:

- To study the characteristics of electron beam when passing through magnetic field and electric field.
- To calculate the ratio between the charge of electron and its mass (e/m).

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
Aim of the experiment:

- Demonstrating the wave nature of electrons through the observation of interference caused by passage of electrons through a polycrystalline graphite lattice (Debye-Scherrer diffraction).
- Determining the wavelength as a function of the anode voltage from the radius of the diffraction rings and the lattice plane spacing of graphite, as well as confirming de Broglie's hypothesis

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
Aim of the experiment:

- NMR experiment is one of classic experiments in recent physics. In this facility one may through the oscilloscope NMR signal of $1H$ and $19F$, measure and calculate g -factor and nuclear magnetic moment value and in turn master and understand NMR principle and phenomena.

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Component
Power supply.
Edge oscillator.
Permeant magnet.
Exploring coil.
Pure water sample.
Water admixed with $CuSO_4$.
PTFE solid sample.
Frequency meter.
Oscilloscope.

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- To observe the phenomenon of electron spin resonance,
- To determine the g -factor of DPPH, and precisely measure the parameters of the geomagnetic field where the system is located, such as the geomagnetic vertical component, geomagnetic horizontal component, and magnetic dip.
- It is built-in microprocessor to automatically determine resonant point and acquire measurement data with high precision.

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
Component
Uniquely designed rotary Helmholtz coil-sensor head combination
Microprocessor-aided auto-determination of resonant point
Computer displaying, processing, and printing of experimental data
Special resonant material for precise measurement
Four-digit display with audio and visual indication of resonant point
Oscilloscope 50 MHz
Digital multi meter UT39A
Cables

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**Aim of the experiment:**

- To observation of Zeeman Effect, enhance the concept understanding of atomic magnetic moment and split quantization in atomic physics.
- To observe spectral lines splitting and their polarization states of 546.1nm mercury atomic spectrum, calculate electron charge-mass ratio based on Zeeman splitting amount.
- To learn the adjustment method to Fabry-Perot etalon.

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Aim of the experiment:

The optical magnetic resonance experiment instrument (shortened as "optical pumping") is used in modern physics experiments.

Involving rich knowledge about physics, such experiments enable students to understand optics, electromagnetism and radio electronics against realistic contexts, and make possible understanding the internal information of atoms qualitatively or quantitatively. They are one of the typical experiments used in spectroscopic teaching. Optical magnetic resonance experiment uses optical pump and photoelectric detecting technology, and thus is a way above ordinary resonance detection technologies in sensitivity. This approach is widely applicable in basic physics research, accurate measurement of magnetic fields, and technical standards making of atomic frequency.

Experiments:

- Observe optical pumping signal
- Measure g-factor
- Measure earth magnetic field (horizontal and vertical components)

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Description	Specifications
Horizontal dc magnetic field	0 ~ 0.2 mt, adjustable, stability $< 5 \times 10^{-3}$
Horizontal modulation magnetic field	0 ~ 0.15 mt (p-p), square wave 10 hz, triangle wave 20 hz
Vertical dc magnetic field	0 ~ 0.07 mt, adjustable, stability $< 5 \times 10^{-3}$
Photodetector	Gain > 100
Rubidium lamp	Lifetime > 10000 hours
High frequency oscillator	55 mhz ~ 65 mhz
Temperature control	~ 90 oc
Interference filter	Central wavelength 795 ± 5 nm
Quarter wave plate	Working wavelength 794.8 nm
Polarizer	Working wavelength 794.8 nm
Rubidium absorption cell	Diameter 52 mm, temperature control 55 oc

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Aim of the experiment:

- The ST365 is a versatile wireless radiation counter engineered to streamline routine sample counting tasks across various fields, including education, training, cutting-edge research, radiation safety, and emerging technologies.
- Its design offers the flexibility to operate as a portable or benchtop device. The ST365 seamlessly interfaces with computers via USB, Ethernet, or WiFi through Spectrum Techniques' STU software. It accommodates both GM and scintillator-type detectors, featuring precise high-voltage adjustment controls.
- The ST365 is offered in two distinct models: Model A, equipped with a 6-digit LED screen, and Model B, boasting a 16-digit backlit LCD.

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Component
ST365 WIRELESS RADIATION COUNTER
GP35 – GM DETECTOR ASSEMBLY
SOFTWARE AND CABLES

***PLEASE NOTE THAT RADIOACTIVE SOURCES ARE NOT INCLUDED**

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Aim of the experiment:

- The unit is a portable training system for studying the principles of solar energy, storage and conversion. Different electrical parameters can be measured due to various digital meters.
- The solar panel and storage battery is connected to a digital charge controller, which is designed with data acquisition and software monitoring function that allows the output from solar panel and battery state to be investigated, as well as controlling the output load.
- Adjustable distance between simulated light source and polycrystal solar cell panel, allowing investigation of system performance inside the laboratory irrespective of outside conditions.

TRAINING CONTENT

Solar cell open circuit voltage test	Solar cell module output characteristics test
Solar cell open circuit current test	Solar cell direct charging experiment
Solar cell maximum power output test	Solar street lamp experiment
Solar cell I-V character test	Solar street lamp charging test
Solar cell P-V character test	Solar street lamp charging/discharging test
Influence of environment to solar cell photovoltaic conversion	Solar fan experiment
Direct resistive load current and voltage test	Solar fan charging test
Direct inductive load current and voltage test	Solar fan charging/discharging test experiment

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CONFIGURATION

Photovoltaic system
Instrument unit
Load unit
Charge controller
Power supply
Inverter grid off
Storage battery
Experiment box and panel:
Accessory

CONFIGURATION SPECS.	
Solar controller	(with electronic manual) with temperature sensor interface and 485 communication interface.
Photovoltaic interface	internally connected to the peripheral solar panel.
Load potentiometer	The resistance value is 0 to 110 ohms, and the power is 2W.
Controller interface	Internally connected to the solar controller
Circuit schematic diagram	an overview of the principle of the system
Brightness adjustment of spotlights	Adjust the brightness of spotlights.
Spotlight start and stop button	2 spotlights can be started independently
Display instrument	There are AC voltmeter 0 to 450V, AC ammeter 0 to 3A, DC voltmeter 0 to 50V, and DC ammeter 0 to 5A.
DC load	here are DC12V fans, DC12V LED lights
DC load light	LED light DC5V
AC load	here are AC220V fans, AC220V lights.
Switch	Set up to turn on/off the AC220V load.
DC motor	The voltage of the DC motor is DC12V
DC step-down module	provide DC5V voltage for the load, USB output charging voltage DC5V
Inverter	convert the input DC12V direct current into AC220V alternating current and provide 100W power to the load.
Power circuit breaker	used to cut off the main power of the control system
Battery	used in conjunction with the solar controller to store and release electrical energy
Spotlight (Simulated light source) mounting plate	Install two spotlights together
6 degree of freedom bracket (Flexible light support)	In order to adjust the orientation of the spotlight, the solar panel can achieve better illumination effect.
Support frame	fix the spotlight bracket, and can also adjust the left and right, front and rear positions properly, and cooperate with the 6-DOF bracket to achieve a better position, which is easy to disassemble and assemble.
Spotlight	The single power of the Spotlight is 500W, and the two are 1000W, which are set up to simulate sunlight.
Solar panels	The output voltage of a single solar panel here is 18V, the current is 0.55A, and the rated power is 10W. There are 2 in total.

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
Aim of the experiment:

- Understand the working principle of solar cells;
- Measure the characteristic parameters such as volt-ampere characteristic curve, open circuit voltage, short circuit current, maximum output power and fill factor of solar cells;
- Understand the working principle of proton exchange membrane electrolysis cell (PEMWE);
- Understand the working principle of proton exchange membrane fuel cell (PEMFC);
- Measure the volt-ampere characteristic curve, open circuit voltage, short circuit current, maximum output power and conversion efficiency of the fuel cell;
- Observe the energy conversion process: light energy → solar cell → electric energy → electrolytic cell → hydrogen energy → fuel cell → electric energy.

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Aim of the experiment:

- Solar power generation
- Demonstrate the process of hydrogen production from water by proton exchange membrane electrolysis
- Demonstrate the process of proton exchange membrane fuel cell power generation
- The results show that the electrochemical reaction of electrolyzed water and power generation is reciprocal

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Component
Solar Module
PEM Electrolyser
Hydrogen And Oxygen Storage Tanks
PEM Fuel Cell
Electric Load; Mounted On A Black Base Plate

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Aim of the experiment:

- The Renewable Energy Science Kit demonstrates the workings of a clean energy technology system on a miniature scale. Power an electrical circuit by solar panel or a wind turbine with profiled blades based on NASA aeronautics. Generate hydrogen through water electrolysis and convert it into electricity using a PEM fuel cell. Whichever combination of technologies you want to explore, this science kit is a comprehensive introduction to the principles behind renewable microgrids.

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Component	
Wind turbine body	Inner Gas containers
Rotor head for profiled blades	Circuit board module base
3 profiled blades for turbine	100-ohm Variable Resistor module
Rotor head adapter for sheet blades	1 Watt Solar panel
3 polypropylene sheet blades for turbine	Adaptors, tubing clincher & purging valve
Rotor unlocking tool	Assembly instructions
Turbine Support base	CD with curriculum manuals
Aluminum wind turbine post	Water/gas tank module base
PEM Electrolyzer	Flexible 2mm banana connecting leads
PEM Electrolyzer base	Transparent silicon tubing
PEM Fuel cell	Plastic plug pins for electrolyzer
PEM Fuel cell base	Battery pack with connecting leads
Hydrogen tank	Syringe
Oxygen tank	

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**Aim of the experiment:**

- A set of equipment permitting quick and easy experiments to demonstrate the effect of greenhouse gases on the absorption of infra-red radiation.

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Component
Base plate
Lamp holder with reflector lamp
Cuvette on stem
Black metal disc
Metal tube, simple
Metal tube, with taps
2 Mounting stems
Silicone hose
Storage case

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Aim of the experiment:

- The unit demonstrates the practical application of such a combined system for student to conduct the investigation of wind power generation.
- Three phase AC permanent magnet generator is used for the wind power generation, which is driven by a DC gear motor to simulate the wind turbine working under different wind speed. By this way to achieve system performance inside the laboratory irrespective of outside conditions.
- Maximum power point tracking type wind controller is adopted for the high efficiency power generation control. Multifunctional meter is adopted for system voltage / current monitoring, with LCD screen.
- Dedicated instrumentation includes AC/DC voltmeter and ampere meters for investigation of system efficiency and characteristics.

TRAINING CONTENT

Wind power generator simulated wind speed experiment	Wind controller resistance load experiment
Wind power generator output voltage experiment	DC resistive and inductive load experiment
Wind controller braking experiment	AC resistive and inductive load experiment
Wind power generator V/I output characteristic curve experiment	Integrated load experiment
Storage battery charging experiment	


Experiment component

Frame structure	- Complete aluminium profile as support frame for wind turbine and simulated wind source - Engraved and colored printed panel, be inserted into grooved aluminium profile - Overall size: 950mm x 455mm x 555mm (approx.)
Wind power generation system	- Wind power generator: three phase AC permanent magnet generator, power 100W, 12VAC, speed 750r/min - Driving motor: DC gear motor 69W, 3000r/min, reduction ratio 1:3, speed adjustable - Speed sensor: 12VDC, NPN type - Wind speed control knob - Wind turbine control buttons and indicators
Instrument unit	- AC voltage meter 0-450V - DC voltage meter 0-50V - AC current meter 0-3A - DC current meter 0-5A - Multifunction meter: LCD display voltage / current / power, with on/off switch
Load unit	- DC12V inductive load: sleeve bearing type fan, 0.2A/4500RPM with protective cover - AC220V inductive load: dual ball bearing type fan, 1.5W/4700RPM with protective cover and ON/FF switch - DC12V resistive load: angle adjustable 5 W LED white lamp - AC220V resistive load: angle adjustable 3W LED white lamp with ON/OFF switch - Variable resistor load: 0Ω...30Ω, 100W adjustable - DC motor 12V: with rotary disc - Pilot lamp DC 12V
Maximum power point tracking wind controller	- Rated voltage 12V, 100W - System related voltage: 12V/24V automatic identification - Charging mode: MPPT - Protection grade: IP67 - Standby current: 3.6mA
Storage battery	12V8AH maintenance free sealed lead acid battery
Accessory	Accessory - Safety whole sealing cable: Conductor cross section: 16AWG tinned copper wire, 3kV, 20A, with axial socket - Power cord

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
**Aim of the experiment:**

- Measure the volt-ampere characteristics of the solar cell when it is forward biased.
- Measure the output characteristics of solar cells under light conditions, as well as open circuit voltage V_{OC} , short circuit current I_{SC} , maximum output power P_{max} , and fill factor FF .
- The light effect of solar cells.

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**Aim of the experiment:**

- An improved model of our plate solar thermal collector, for both qualitative and quantitative experiments. The solar thermal collector is supplied with a circulation pump and hot water container. The hot water container can also be used empty to illustrate the heating of air by a radiator. Fitted with a peristaltic pump that easily fills and empties pipes and hoses.

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Aim of the experiment:



The X-ray apparatus allows you to undertake a wide variety of experiments on the following topics:

- Properties of X-rays: (Penetration by X-rays / straight-line propagation / Ionization / X-ray photography)
- Dosimetry and radiation protection
- Fluorescence Refraction of X-ray beams: (Laue imaging / Debye-Scherer imaging)
- Shielding from X-rays
- Bragg reflections
- Absorption experiments
- Duane-Hunt displacement law (determination of Planck's constant h)
- Distance law
- Moseley's law

Experiment content:

	Component
1	Fine adjustment trimmer for goniometer
2	Trimming potentiometer for setting emitter current.
3	High voltage ON/OFF button
4	Housing
5	Base
6	Timer
7	Key switch for mains supply
8	Jack socket for emitter current measurement, giving an output voltage proportional to current
9	Cover
10	Screening plate with radiation warning label
11	Dumbbell shaped slot for locking the cover
12	Protractor scale for measuring arm
13	Measuring arm with magazine for photographic slides
14	Pressure plate
15	Protractor scale for sample holder
16	Sample holder
17	Radiation exit window with brass collimator
18	Securing screws for leaded glass dome
19	X-ray tube
20	Access holes
21	Leaded glass dome
22	Locking pin for cover
23	Length scale for focal point
24	4-mm sockets for addition of driver motor
25	High voltage selector switch
26	Mains voltage indicator lamp (white)
27	High voltage indicator lamp (red)
28	Circuit breaker for mains supply
29	Circuit breaker for high voltage supply
30	Mains voltage selector switch
31	Ventilation grille
32	Mains socket

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Aim of the experiment:

- This set helps to understand the Doppler Effect and its application in medicine for blood flow measurements. With the ultrasonic pulse Doppler „Flow Dop” and a special Doppler probe ultrasonic waves can be generated, and frequency shifted scattering signals can be measured. With the Doppler prisms various angles of incidence can be realized. In this way the relation between frequency shift of the Doppler Effect, the angle of incidence, transmission frequency and flow velocity can be determined. With a flow profile the occurrence of a laminar or turbulent flow can be measured. A realistic arm model is used to simulate the application of the Doppler Effect in medicine. With a Doppler sonography the influence of a stenosis on the flow profile can be investigated. A pump generates different flow types (continuous and pulsatile) and can simulate the human blood circulation. The measured Doppler signals can be presented acoustically as well as in a color-coded Doppler spectrum, whereas the results and images are not much different to measurements of those on patients.

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Aim of the experiment:

- Master the measurement method of hearing and hearing Threshold.
- Determine the hearing threshold curve of the human ear.

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Description	Specifications
Signal source	Frequency range: 20 ~ 20 khz. standard sine wave (smart key controlled)
Digital frequency meter	20 ~ 20 khz, resolution 1 hz
Digital sound strength meter (db meter)	Relative -35 db to 30 db
Headset	Monitoring grade
Power consumption	< 50 w
Instruction manual	Electronic version

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Aim of the experiment:

- Understand the working principle of the gas pressure sensor and test its characteristics.
- Use gas pressure sensor, amplifier and digital voltmeter to construct a digital pressure gauge and calibrate it with a standard pointer pressure gauge.
- Understand the principle of measuring human heart rate and blood pressure, use pulse sensor to measure pulse waveform and heartbeat frequency, and use the constructed digital pressure gauge to measure human blood pressure.
- Verify boyle's law of the ideal gas. (optional)
- Use slow scanning long afterglow oscilloscope (need to be purchased separately) to observe the body pulse waveform and analyze the heart beat, estimate heart rate, blood pressure and other parameters. (optional)

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Aim of the experiment:

- Study the braking reaction time of cyclist or car driver when the signal light is changed.
- Study the braking reaction time of cyclist when hearing the sound of a car horn.

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Description	Specifications
Car horn	Volume continuously adjustable
Signal light	Two sets of led arrays, red and green colors respectively
Timing	Accuracy 1 ms
Time range for measurement	Unit in second, signal may appear randomly within the set time range
Display	LCD module

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Aim of the experiment:

- The Abbe Refractometer is intended for measuring the refractive index nD and the average dispersion nF – nC of transparent or translucent liquid, or solid substance (of which mainly intended for measuring transparent liquid). In case of the instrument being connected with a thermostat, it will be able to measure the refractive index nD at a temperature between 0° and 70°.

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
Head	Monocular Head
Refractive index ND measurement range	1.300~1.700
Refractive index ND measurement accuracy	0.0002
Refractive index ND Min. division	0.001
Glucose concentration (%) measurement range	0~95
Glucose concentration (%) Min. division	0.5

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
Aim of the experiment:


- Measure Muon Lifetime
- Demonstrate Relativistic Time Dilation
- Measure Local Muon Flux
- Measure Sea Level Muon Charge Ratio
- Convenient Source of Genuinely Random Numbers
- Create Simulated "Muons" and measure their Lifetime
- Study Processing of Photomultiplier Signal

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Component	Specification
Three "modules" comprise the Muon Physics hardware:	
The detector	Detector Size: Diam. 16.5 cm Height 35.5 cm.
The electronics	PMT: 10-stage bialkali photocathode Diam. 5.1 cm.
User supplied computer.	Timing FPGA: Bin Size (resolution) 20 ns.
	Dynamic Range 20 μ sec.
	Timing clock frequency: 50 MHz \pm 5 kHz.
	Power Consumption (excluding PC): 25 watts.
	Typical detected muon decay rate: 1 event/minute.
	Supported I/O port protocols: Serial and USB.

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Aim of the experiment:

- This experiment adopts the principle of sound source positioning in acoustic emission technology, simulates satellites with ultrasonic signal transmitters, and simulates user GPS receivers with ultrasonic signal receivers.
- Experimental simulation of two-dimensional plane positioning and intuitive three-dimensional space GPS positioning.


OTHER EXPERIMENTS

Electrochemical
Equivalent of Copper

Electromagnetic
Induction
Demonstration


Air Flow Apparatus


Centripetal Force
Apparatus


Hoffman Electrolysis
Apparatus

Hydroelectric Power
Plant As A Model

Internal Resistance
Of Ammeter &
Voltmeter


Jet Nozzle (Flow
laminator)


Torsion Balance

Rigidity Modulus
Experiment

Trans-Ratio
Experiment

Demonstration
Transformer

Torsion Balance

Light Energy
Apparatus

Velocity of Sound
Using Tuning Fork

Na lamp (Sodium Lamp)

A low-pressure sodium lamp can emit a strong 589.0nm and 589.6nm, which is a good monochromatic and is often used as a monochromatic light source in polarizers, refractometers, polarimeters, etc.

SPECIFICATION OF THE BULB LAMP:

Power (w): 20W

Average working life: 1000hr

Bulb diameter (mm): 28

Length (mm): 155

Illumination center (mm): 75

Hg lamp (Mercury Lamp)

A low-pressure mercury lamp can emit a strong 253.7 nm, 404.7nm, 435.8nm, 546.1nm, 577nm, 579nm of mercury lines, often used as a monochromatic light source in optical rotator, refractometer, polarimeter and other instruments.

SPECIFICATION OF THE BULB LAMP:

Power (w): 20W

Average working life: 1000hr

Bulb diameter (mm): 28

Length (mm): 155

Illumination center (mm): 75

D2 Lamp (Deuterium Lamp)

UV spectral lamp with long life (2000 hours on average), high stability and high intensity. This lamp is used for various analytical instruments as UV Spectro light source. This lamp can be used for: Spectrophotometer, High performance liquid chromatography, atomic absorption spectrophotometer. (background correction) and Film thickness meter Output variation (30 minutes after lamp ignition):

Drift (variation of light output with passage of time): $\pm 0.5\%/h$

Fluctuation (Max variation of light output in short period): 0.1% P-P

Average life: 2000 hours

Brand: OSRAM

Cd-10 Lamp (Cadmium Lamp)

UV spectral lamp with long life (2000 hours on average), high stability and high intensity. This lamp is used for various analytical instruments as UV Spectro light source. This lamp can be used for: Spectrophotometer, High performance liquid chromatography, atomic absorption spectrophotometer. (background correction) and Film thickness meter Output variation (30 minutes after lamp ignition):

Drift (variation of light output with passage of time): $\pm 0.5\%/h$

Fluctuation (Max variation of light output in short period): 0.1% P-P

Average life: 2000 hours

Brand: OSRAM

SPECTRUM TUBES WITH POWER SUPPLY

SPECTRUM TUBES

Spectrum Tubes demonstrate the wave properties of different gases.

- Ideal for Wavelength Determination
- For Introducing Students to Better-Known Gases
- Glass Tube Length: 10.25" (26cm)
- Capillary Portion Length: 4" (10cm)
- Tubes Work in a Variety of Spectroscopic Applications

Ideal for wavelength determination and introducing students to the spectra of better-known gases, these glass tubes are 10.25" (26cm) long, and the capillary portion is about 4" (10cm) long. The spectrum tube assemblies can be used in a variety of spectroscopic applications.

When they are observed through a spectroscope, the emission lines corresponding to the energy transition levels of the gaseous atoms or molecules appear as distinct components.

NOTE: These tubes are designed to operate 30 seconds ON - 30 seconds OFF.

SPECTRUM TUBES LIST WE OFFER:

- Spectrum Tube Ar (Argon)
- Spectrum Tube H₂ (Hydrogen)
- Spectrum Tube He (Helium)
- Spectrum Tube Hg (Mercury)
- Spectrum Tube Kr (Krypton)
- Spectrum Tube Ne (Neon)

SPECTRUM TUBE POWER SUPPLY

A specially designed box holds the tube firmly and prevents the student from touching the electrodes. A black panel behind the Spectrum tube eliminates distracting ambient light and protects the tube from breakage. Operates at 220V A/C at 50Hz and is suitable for 20/26 cm tubes.

Also, the power supply can be used with a timer to maintain the lamp and extend its operating lifetime as These tubes are designed to operate 30 seconds ON - 30 seconds OFF.

TECHNICAL DETAILS:

- Input Voltage = 220VAC, 50Hz
- Output Rated: 5,000V without spectrum tube 1,000V @10mA with spectrum tube.
- Operating Temperature: 10°C to 40°C

LASER SOURCES

HELIUM-NEON LASER SOURCES

This He-Ne laser is reasonable and fashionable designed, characterized by stable and reliable performance, good mono-chromaticity and directivity, high brightness, long life, good adapted to electric network and widely applied in many fields.

SPECIFICATIONS:

Wavelength: 632.8nm
Output laser power: >2mW.
Continuously working time: >8hours
Dimensions: 380*90*100mm

Minimum output power: 1.6mW
Beam diameter (TEM₀₀): 0.58mm
Beam divergence (TEM₀₀): 1.4mrad
Polarization: Random
Output stability: 5%.
Wavelength: 632.8nm
Transverse mode: TEM₀₀ or TEM₀₁

DIOD LASER SOURCE

The laser module emits a light in the form of a red light. The maximum output power is < 1 mW. Voltage supply from 3 to 5 V

SPECIFICATIONS:

Wavelength: 650 nm
Maximum output power: < 1 mW.
Supply current: 40 mA.
The knob to adjust the light intensity (the focus of the laser dot).
Supply voltage: 5V2A.
Size: 120*80*40mm

LASER POINTER

Easy to use Pen diode laser pointer. Runs on a 3.7v chargeable lithium battery. This laser is monochromatic and has a very well-defined wavelength. Can be used in interference experiments as well as in geometrical optics. Optical power is about 1mW.

LASER COLORS:

- LASER Pointer RED
- LASER Pointer BLUE
- LASER Pointer GREEN

LASER RAY BOX

This laser ray box is an excellent classroom tool for demonstrating optical ray paths to students. Teaching and studying geometric optics have never been easier with this ray box. Just add a few lenses and mirrors, and you are ready to go. This model comes with five red lasers, which you can control with a provided switch to turn either one, three, or all five beams at a time. It also comes with a 6V DC power adapter, so you don't have to use batteries like other ray boxes.

SPECIFICATIONS:

- Laser Output Power: <1 mW Class IIIa (each laser beam)
- Wavelength: 650 nm
- Spacing between lasers: 13 mm
- Toggle between 1, 3, or 5 beams.
- On/Off switch
- Ray box size: 85 mm x 50 mm x 24 mm

POWER SUPPLIES

LINEAR DC POWER SUPPLY 0 - 30VDC 0-5A

series single-channel linear DC power supplies are economical benchtop tools for DIY electronics hobbyists, electronics manufacturing and testing, and other applications.

SPECIFICATIONS:

- Output voltage: 0~30V.
- Output current: 0~5A.
- Output power: 150W
- Load regulation: $CV \leq 0.01\% + 5mV$, $CC \leq 0.1\% + 10mA$.

ARABIA VARIABLE (0 - 20) VDC (1A) POWER SUPPLY.

The Arabia DC variable power supply is ideal for use as an economical classroom or laboratory power supply and will provide years of trouble-free service.

SPECIFICATIONS:

- Output Ports (+ / -): 4mm.
- Digital panel display (amperes and voltage)
- Dimensions: 25cm*19cm*9cm.
- Useful as a Battery Eliminator
- Output Voltage: 0 - 15V DC @ 1A max.
- Output Voltage Selection: Variable (Coarse) DC.
- Output Power: 15W.
- Supply Voltage: AC 220V \ 50Hz.

ARABIA AC SELECTED POWER SUPPLY (3-6-9-12-15) V 3A

ARABIA DC SELECTED POWER SUPPLY (4-8-12-16) V 3A

The Arabia DC/AC selected power supply is ideal for use as an economical classroom or laboratory power supply and will provide years of trouble-free service.

SPECIFICATIONS:

- Output Ports (different values / GND): 4mm.
- Dimensions: 16cm*16cm*13.5cm.
- Output Voltage: (3-6-9-12-15) V AC @ 3A max.
- Output Voltage: (4-8-12-16) V DC @ 3A max.
- Output Voltage Selection: Variable (selector) DC\AC.
- Output Power: 9W-50W.
- Supply Voltage: AC 220V \ 50Hz.

POWER SUPPLY ADAPTOR

SPECIFICATIONS:

- Output: DC 12V \ 5A
- Input: 240VAC, 50 1.8A
- Out diameter of the connector: 5.5mm
- Inner diameter of the connector: 2.1mm

ENVIRONMENTAL MEASUREMENT TOOLS

DIGITAL ANEMOMETER

stable, safe and reliable split-type anemometer, featuring quick detection, low consumption and high accuracy. It can be widely used in mining, electric, iron and steel, petrochemical, energy-saving, navigation, fan manufacturing, exhaust ventilation, sporting and more sectors.

Specifications	Range	Resolution	Accuracy
Wind speed	0.4~30m/s	0~99: 0.01 100~999: 0.1 ≥1000: 1	±(5%+0.5m/s)
	1.4~108km/h		±(5%+1.5km/h)
	78~5905ft/min		±(5%+100ft/min)
	0.7~58knots		±(5%+1knots)
	0.8~67mph		±(5%+1mph)
Temperature	-10~5°C	0.1°C	±2.0
	14~122°F	0.2°F	±4.0
Wind scale	0~12	1	±1

DIGITAL LIGHT METER

a stable, safe, reliable mini digital light meter, widely used in light manufacturing, photo studios, interior renovation, construction, greenhouses.

Specifications	Range	Resolution	Accuracy
Illuminance	0~199,900Lux	1Lux (0~9999Lux)	±(4%+8)
		10Lux (≥10,000Lux)	±(5%+10)
		100Lux (≥100,000Lux)	±(5%+10)
		0.1Fc (0~999.9Fc)	FC=Lux/10.76
		1Fc (≥1000Fc)	
		10Fc (≥10,000Fc)	

DIGITAL TEMPERATURE HUMIDITY METER

a stable, safe and reliable mini digital temperature humidity meter. With split type design, it is convenient and fast to measure the temperature accurately in different environments. The split type design is perfect for monitoring hard to reach and small ventilation ducts.

Specifications	Range	Resolution	Accuracy
Temperature (°C)	-10°C~60°C	0.1°C	±1.0°C
Temperature (°F)	14°F~140°F	0.2°F	±2.0°F
Humidity	0%RH~100%RH	0.1%RH	±5% (20~80%RH)

SOUND LEVEL METER

A sound level meter can measure noise level and output results in db. A or C weighted modes can be selected depending on the application. UT352 consumes very little power and can continuously work for up to 20 hours. They are suitable for industrial and environmental sound level measurement.

Specifications	Range	UT352
Range	30~80dB	±1.5dB
	50~100dB	±1.5dB
	60~110dB	±1.5dB
	80~130dB	±1.5dB
Frequency	31.5Hz~8000Hz	√
Frequency Weighting	A&C	√
Sample Rate	FAST: 8 Times/s	√
	SLOW: 1 Times/s	√

OPTICAL COMPONENTS

LENS, DOUBLE CONVEX: Explore our selection of precision-engineered biconvex lenses, designed to facilitate a wide range of geometric optics experiments with accuracy and clarity. These lenses are indispensable tools for educators, researchers, and enthusiasts seeking to explore the fundamental principles of light and optics.

DIAM, MM	FOCAL LENGTH, MM	POWER
50	50	+20D
50	100	+10D
50	150	+6.6D
50	200	+5D
50	250	+4D
50	300	-

LENS, DOUBLE CONVEX: Explore our selection of precision-engineered biconvex lenses, designed to facilitate a wide range of geometric optics experiments with accuracy and clarity. These lenses are indispensable tools for educators, researchers, and enthusiasts seeking to explore the fundamental principles of light and optics.

DIAM, MM	FOCAL LENGTH, MM	POWER
50	100	-10D
50	200	-5D
50	250	-4D

CONCAVE AND CONVEX MIRROR

Optically Worked, back silvered, with a protective coating

DIAM, MM	FOCAL LENGTH, MM
50	100
50	300

CROWN PRISM: The material of this product is glass. It is made of optically worked crown glass which has a relatively low refractive index of approximately 1.52 and a low dispersion. The crown glass prisms are mainly used with spectrometers. Students can also understand how chromatic aberrations can be reduced by combining a flint glass prism with a crown glass prism.

Available with two sizes


- 32*32mm
- 25*25mm

FRESNEL BIPRISM: The material of this product is glass. It is made of optically worked crown glass which has a relatively low refractive index of approximately 1.52 and a low dispersion. The crown glass prisms are mainly used with spectrometers. Students can also understand how chromatic aberrations can be reduced by combining a flint glass prism with a crown glass prism.

Available with two sizes

- 32*32mm
- 25*25mm


GLASS PRISM AND RIGHT-ANGLE PRISM: A prism is a three-dimensional solid object with two identical ends, flat faces, and equal cross-sections.

- 50mm*50m*50mm

EDUCATIONAL DIFFRACTION GRATING: A Diffraction Slide for illustrating the principles of diffraction grating and how its properties relate to the number of lines per unit length. These gratings have line densities of 100 lines/mm, 300 lines/mm, and 600 lines/mm, respectively, and are safeguarded by protective cover glasses

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